MASTERMIND® Q.F.T.®

- ITA** Sistema integrato per l'esecuzione di test di filtrabilità e test SDI
- ENG** Integrated system for performing filterability and SDI tests
- DEU** Integriertes System zur Durchführung von Filtrierbarkeits- und SDI-Tests
- FRA** Système intégré pour la réalisation d'essais de filtrabilité et SDI
- ESP** Sistema integrado para realizar pruebas de filtrabilidad y prueba SDI

38cm x 48cm (14.96 in x 18.90 in)

Misura
Size
Größe
Taille
Medidas

MASTERMIND® Q.F.T.®

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MASTERMIND® Q.F.T.®

Integrated system for performing filterability and SDI tests.

DESCRIPTION

Mastermind® Q.F.T.® is a laboratory instrument capable of **performing both filterability indices on wine and fouling indices on water**: by means of an accurate weighing test, it calculates indexes on wine and water in complete automatism, with temperature control.

Mastermind® Q.F.T.® calculates these filterability indices automatically, based on temperature-controlled weight readings.

The results, which can be repeated several times due to the high degree of standardisation allowed by the instrument, are displayed on a PC display (not supplied) and stored in a database (MS Access), or printed in a report (or stored in pdf format).

A dedicated software is included. The heart of the instrument is a precision scale developed specifically for this application in collaboration with Gibertini Elettronica Srl, a world leader in the production of precision weighing instruments.

Mastermind® Q.F.T.® also records the temperatures at which the test is carried out: this allows a true comparison of the data collected by the instrument at even very different times. It is also possible to modify the test parameters to perform filterability indices in special applications (e.g. juices). The instrument performs the washing through its own work cycle.

STRENGTHS

- High accuracy and repeatability
- Intuitive management software
- Measurement of IF, MFI and Vmax parameters in a single test
- Highly accurate measurement based on weight reading
- Possibility of checking the clogging effect of service waters (S.D.I. test)

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JUCLAS S.r.l.