MASTERMIND® REMOVE LAB UNIT

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MASTERMIND® REMOVE LAB UNIT

Laboratory dealcoholisation unit.

- ⊕ Alcohol reduction without altering the colour and aromas
- ⊕ Faithful reproduction of the industrial process
- ⊕ Delicate process at room temperature and pressure

DESCRIPTION

The **MASTERMIND® REMOVE LAB UNIT** dedicated to alcohol reduction is designed to reproduce in a faithful and controlled manner the process used in the industrial plants of the **MASTERMIND® REMOVE** line. It is a tool designed for research activities, preliminary tests, protocol development and qualitative assessments, which allows small volumes of wine or other alcoholic beverages to be processed prior to large-scale production.

The operating principle is based on **perstraction**, the same patented technology used in industrial plants (patent no. 0001387862). The wine passes through a **contactor membrane with a hydrophobic interface**, which allows the selective passage of ethanol in gaseous form to the **extraction solvent** (water). The process takes place at **room temperature and pressure**, without oxygen enrichment and without changing the colour, structure, or aromatic complement of the wine. The aim is to achieve a controlled reduction in alcohol content while preserving the initial organoleptic profile, so as to explore the wine's expressive 'sweet spot' experimentally.

From an operational point of view, the laboratory version allows the **extraction of up to 1 L/h of alcohol**, which is adequate for testing requirements on small product volumes. The unit is capable of handling **volumes between 3 and 30 L** and is designed to provide instructive and repeatable experimental evaluations; the alcohol reduction achievable in one cycle typically ranges from **0.1% vol to total dealcoholisation**, depending on the duration of the treatment. The operation is simple and completely manual, ensuring the operator has maximum control over flow rates, timing, and process evolution.

MASTERMIND® REMOVE LAB UNIT is an agile, compact, and scientifically reliable tool for the study of alcohol reduction. It allows the formulation of protocols that can be replicated on larger capacity industrial plants.

The use of the device requires minimal amounts of energy and water, making it ideal for both preliminary research and pilot sampling for the development of low or zero-alcohol products, an area that has also become increasingly relevant at a regulatory level.

STRENGTHS

- **Compact, intuitive, and simple-to-use** instrument
- **Faithful reproduction** of the MASTERMIND® REMOVE industrial process
- Process conducted at ambient temperature and pressure
- **High scientific reliability** for repeatable tests
- No connection to the electrical grid required

AREAS OF APPLICATION

- **Research and development** in the field of wine dealcoholisation
- Definition of the alcoholic '**sweet spot**' through controlled trials
- Study of organoleptic profiles with different alcohol contents
- **Pilot trials** and preliminary tests prior to large-scale production
- Quality assessments for wines intended for partial or total alcohol reduction

MANUAL MODELS

- MMR LAB Unit — up to 1 L/h of extracted alcohol

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JUCLAS S.r.l.