SMARTCHECK®

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SMARTCHECK®

Technology 4.0 for assessing the tartaric stability of wines.

- + Monitoring of stabilisation processes.
- + Cellar process decision guide.
- + Improved ability to measure and monitor production inputs.

DESCRIPTION

SMARTCHECK® is an innovative instrument designed and manufactured thanks to the interdisciplinary expertise within the **VASONGROUP** research laboratories. The instrument is capable of performing an evaluation of the tartaric stability of wines by means of an analysis based on the principle of the **mini-contact test**, which involves following the fall in conductivity (crystal growth) following the addition of potassium bitartrate to a wine subjected to continuous agitation at low temperature. The analysis with **SMARTCHECK®** can be conducted at 5 different temperatures for the time required by the operator; this makes it a **"direct" instrument** that can be used to carry out more in-depth scientific studies and to run the production stages, by obtaining a specific parameter, even in terms of repeatability. **The results can be viewed in real time on an intuitive three-screen display**, where you can follow the test via temperature indicators and initial, instantaneous and differential conductivity parameters.

STRENGTHS

- The test is performed with a choice of 5 possible temperatures.
- The test can be performed for any length of time.
- Small in size, to be easily placed on a laboratory counter.
- Electronics area protected from condensation and well separated from the fluid area.
- Easy washing of the fluid management area.
- Modulated and automatic precision cooling to 0.1 °C (32.18 °F), in a water stream.
- Setup with 4.0 wine cellar management system.
- Original JUCLAS software.

AREAS OF APPLICATION

- Parameterisation and analysis for each wine.
- A guide for decisions regarding winemaking processes.
- Better measuring and monitoring capacity of production input.
- Improvement of internal efficiency.
- Monitoring of stabilisation processes and subsequent optimisation.
- Production of characterising databases for the various batches and vintages.

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JUCLAS S.r.l.